

Anti-Carboxysome shell protein CcmK4 antibody

Catalog: PHY5334S

Product Information

Description: Rabbit polyclonal antibody

Background: ccmK4
Synonyms: ccmK4

Immunogen: KLH-conjugated synthetic peptide (18 aa from C terminal section) derived from

Synechocystis sp. PCC 6803 ccmK4 (slr1839).

Form: Lyophilized

Quantity:150 μgPurification:Serum

Peptide affinity form antibody available upon request at info@phytoab.com.

Reconstitution: Reconstitution with 150 µl of sterile water.

"Note: please spin tube briefly prior to opening it to avoid any losses that might

occur from lyophilized material adhering to the cap or sides of the tube".

Stability &Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

Storage: 12 months from date of receipt, -20 to -70°C as supplied.

6 months, -20 to -70°C under sterile conditions after reconstitution.

1 month, 2 to 8°C under sterile conditions after reconstitution.

Shipping: The product is shipped at 4°C. Upon receipt, store it immediately at the

temperature recommended above.

Application Information

Recommended Dilution: Western Blot (1:1000-1:2000)

Note: Optimal dilutions/concentrations should be determined by the

end user.

Expected / apparent MW: 12 / 11 kDa


Confirmed Reactivity: Synechocystis sp. PCC 6803

Predicted Reactivity: For homologues with other species especially algae, please contact

tech support at tech@phytoab.com.



Application Example



PCC6803+G: 30 µg whole-cell lysate protein from WT of

Synechocystis sp. PCC 6803 (Dark Plus Glucose).

PCC 6803+L/G: 30 µg whole-cell lysate protein from WT of

Synechocystis sp. PCC 6803 (Light Plus Glucose).

PCC 6803+L: 30 µg whole-cell lysate protein from WT of

Synechocystis sp. PCC 6803 (Light without Glucose).

Electrophoresis: (5-20)% Gradient gel

Transfer: blotting to NC (nitrocellulose) membrane for 1 h.

Blocking: 5% skim milk at RT or 4° C for 1 h.

Primary antibody: 1:1000 dilution overnight at 4°C.

Secondary antibody: 1:10000 dilution using Goat Anti-Rabbit IgG

H&L (HRP) (Cat# PHY6000).

Detection: using chemiluminescence substrate and image were

captured with CCD camera.