

Anti-PHOSPHATE TRANSPORTER 4;2 antibody

Catalog: PHY0754S

Product Information

Description:	Rabbit polyclonal antibody	
Background:	PHT4;2 (AT2G38060) is a root plastid-localized Pi transporter and contributes	
	to Pi transport in root plastids.	
Synonyms:	PHT4;2, PHOSPHATE TRANSPORTER 4;2	
Immunogen:	KLH-conjugated synthetic peptide (13 aa from N terminal section) derived from	
	Arabidopsis thaliana PHT4;2 (AT2G38060).	
Form:	Lyophilized	
Quantity:	150 µg	
Purification:	Serum	
	Peptide affinity form antibody available upon request at <u>info@phytoab.com</u> .	
Reconstitution:	Reconstitution with 150 µl of sterile water.	
	"Note: please spin tube briefly prior to opening it to avoid any losses that might	
	occur from lyophilized material adhering to the cap or sides of the tube".	
Stability &	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.	
Storage:	12 months from date of receipt, -20 to -70 $^\circ \!\! \mathbb{C}$ as supplied.	
	6 months, -20 to -70 $^\circ\!\!\!\!\!^\circ$ under sterile conditions after reconstitution.	
	1 month, 2 to 8 $^\circ\!{ m C}$ under sterile conditions after reconstitution.	
Shipping:	The product is shipped at $4^\circ\! \mathbb{C}$. Upon receipt, store it immediately at the	
	temperature recommended above.	


Application Information

Recommended Dilution:	Western Blot (1:1000-1:2000)
	Note: Optimal dilutions/concentrations should be determined by the
	end user.
Expected / apparent MW:	55 / 50 kDa
Confirmed Reactivity:	Arabidopsis thaliana
Predicted Reactivity:	For more species homologues information, please contact tech
	support at <u>tech@phytoab.com</u> .

Research Use Only



Application Example



TP: 20 µg, 40 µg and 80 µg total protein from *Arabidopsis thaliana*, respectively.

Electrophoresis: 12% SDS-PAGE

Transfer: blotting to NC (nitrocellulose) membrane for 1 h.

Blocking: 5% skim milk at RT or 4°C for 1 h.

Primary antibody: 1:1000 dilution overnight at 4°C.

Secondary antibody: 1:10000 dilution using Goat Anti-Rabbit IgG H&L(HRP) (Cat# PHY6000).

Detection: using chemiluminescence substrate and image were captured with CCD camera.



Research Use Only

PhytoAB Inc.

Tel:1-800-380-9586 info@phytoab.com www.phytoab.com